




HA1(B/Texas/02/2013)

Product Name	HA1(B/Texas/02/2013)
Catalog No	IT-003-B20p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(B/Texas/02/2013) protein (a.a.16-363) (GenBank #AGO05755).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Note:

E: HA1(B/Texas/02/2013) Protein