




HA1(A/Guangdong/17SF003/2016)(H7N9)

Product Name	HA1(A/Guangdong/17SF003/2016)(H7N9)
Catalog No	IT-003-0075p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(A/Guangdong/17SF003/2016)(H7N9) protein (a.a.19-344) (see below).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ:

DKICLGHHAVSNGTKVNTLTERGVEVVNATETVERTNTPRICSKGKRTVDLGGCGLLGTITGPPQC
DQFLEFSADLIIRREGSDVCYPGKVFVNEEALRQILRESGGIDKEPMGFTYNGIRTNGVTSACRRSG
SSFYAEMKWLLSNTDNAAFPQMTKSYKNTKESPAIIVWGIHHSVSTAEQTKLYGSGNKLVTVGSS
NYQQSFVPSGARPVNGQSGRIDFHWLILNPNDTVTFNFNGAFIAPDRASFLRGKSMGIQSRVQV
DANCEGDCYHSGGTIISNLPFQNIDSRAVGKCPRYVKQRSLLLATGMKNVPEVPRKRRTARG