




# HA1 (A/Hubei/29578/2016)(H5N6)

Product Name	HA1 (A/Hubei/29578/2016)(H5N6)
Catalog No	IT-003-00538p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1 (A/Hubei/29578/2016)(H5N6) protein (a.a.17-339) (see below).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Note:

E: HA1 (A/Hubei/29578/2016)(H5N6) Protein

## Protein SEQ:

DQICIGYHANNSTEQVDTIMEKNVTVTTHAQDILEKTHNGKLCDLNGVKPLILKDCSVAGWLLGNP  
MCDEFIRVPEWSYIVERANPANDLCYPGNLNDYEELKHLISRINHFECTLIIPKSSWPNHTSSGVSA  
ACPYLGKPSFFRNVVWLTKKNDAYPTIKMSYNNTNREDLLILWGIHHSNNAEEQTNLYKNPTTYV  
SVGTSTLNQRVVPKIATRSQVNGQSGRMDFFWTILKPDDAIHFESNGNFIAPPEYAYKIVKKG DSTI  
MKSEMEYGNCNTKCQTPIGAINSSMPFHNIHPLTIGECPKYVKS NKLVLATGLRNSPLRE