


gp140(CM235/CRF01_AE)(HIV-1)



Product Name	gp140(CM235/CRF01_AE)(HIV-1)
Catalog No	IT-001-150ΔTMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	gp140(CM235/CRF01_AE)(HIV-1) (aa 33-657; GenBank# AAG28611) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp140(CM235/CRF01_AE)(HIV-1)

Note: