

# Envelope glycoprotein p57 (Borna disease virus /strain V)



Product Name	Envelope glycoprotein p57 (Borna disease virus/strain V)
Catalog No	IT-028-001p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Envelope glycoprotein p57 (Borna disease virus/strain V) (aa 23-467; GenBank# P52638) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥93%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Envelope glycoprotein p57 (BDV-V)

Note: