




Nucleoprotein (Machupo mammarenavirus)

Product Name	Nucleoprotein (Machupo mammarenavirus)
Catalog No	IT-008-012Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M urea.
Protein	Nucleoprotein (Machupo mammarenavirus) (aa 1-564; GenBank# P26578) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Nucleoprotein (Machupo mammarenavirus)

Note: