


S1 (CoVID-19/Wuhan)



Product Name	S1 (CoVID-19/Wuhan)
Catalog No	IT-002-030p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	The new coronavirus S1 (CoVID-19/Wuhan) (aa 23-692; GIAID# EPI_ISL_402121) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: S1 (CoVID-19/Wuhan)

Note: