



P2C (HAV) (aa 1121-1234)

Product Name	P2C (HAV) (aa 1121-1234)
Catalog No	IT-004-051Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea.
Protein	GST-fused P2C (HAV/Hu/Los Angelos/LA/1975) immune-dominant regions, amino acids 1121-1234 (GenBank Accession No. P06441).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>92% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

