


HA1(A/swine/Jiangsu/J004/2018) (H1N1)



Product Name	HA1(A/swine/Jiangsu/J004/2018)(H1N1)
Catalog No	IT-003-SW17p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA1 (A/swine/Jiangsu/J004/2018) (H1N1)(aa 18-344; GenBank# QEM24690), a G4 EA strain.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazaIT-300-177prdous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(A/swine/Jiangsu/J004/2018)
(H1N1)

Note: