

# gp120(TV1)(HIV-1/Clade C)



---

Product Name	gp120(TV1)(HIV-1/Clade C)
Catalog No	IT-001-170p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp120 (TV1)(HIV-1/Clade C)(a.a. 31-517) (Genebank #: ACF75938)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

