



gp120(BG505.W6M.ENV.C1)(Clade A)

| | |
|----------------|--|
| Product Name | gp120(BG505.W6M.ENV.C1)(Clade A) |
| Catalog No | IT-001-176p |
| Source | Viral protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | 6xHis tagged HIV-1 gp120(BG505.W6M.ENV.C1)(Clade A) protein (a.a.30-508) (Genbank No. ABA61515) |
| Applications | WB standard, antibody ELISA, etc |
| Purity | ≥ 90% |
| Storage | Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

