




# gp120(ZM197M.PB7)(HIV-1)

---

Product Name	gp120(ZM197M.PB7) (Clade C)
Catalog No	IT-001-RC4p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged HIV-1 gp120 (ZM197M.PB7)(Clade C) protein (a.a.33-525) (Genebank No. DQ388515).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(ZM197M.PB7) protein

Note: