




gp120(ZM214M.PL15)(HIV-1)

Product Name	gp120(ZM214M.PL15) (Clade C)
Catalog No	IT-001-RC5p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (ZM214M.PL15)(Clade C) protein (a.a.40-510) (Genebank No. DQ388516).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(ZM214M.PL15) protein

Note: