




HA₁ (H2N2)

Product Name	HA ₁ (H2N2)(A/Canada/720/05 (H2N2))
Catalog No	IT-003-003p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged HA ₁ (H2N2)(A/Canada/720/05 (H2N2)) protein (a.a.13-340) (Genebank accession #AAZ28987)
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified HA₁ (H2N2) protein