




## HA1(H3N2)(A/Philippines/2/82)

---

Product Name	HA1(H3N2)(A/Philippines/2/82)
Catalog No	IT-003-00416p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(H3N2)(A/Philippines/2/82)(17-345aa; GenBank# ADJ41805)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(H3N2)(A/Philippines/2/82)

Note: