


HA1(A/Hawaii/22/2012)(H3N2)



Product Name	HA1(A/Hawaii/22/2012)(H3N2)
Catalog No	IT-003-00426p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(A/Hawaii/22/2012)(H3N2) (a.a.17-345) (same as HA1(A/Quebec/019/2012(H3N2)), Genbank #: AGD99012).
Applications	Protein standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(A/Hawaii/22/2012)(H3N2)