



HA1(A/Switzerland/9715293/2013)(H3N2)

Product Name	HA1(A/Switzerland/9715293/2013)(H3N2)
Catalog No	IT-003-00428p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(A/Switzerland/9715293/2013)(H3N2)(see below for protein sequence)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(A/Switzerland/9715293/2013)(H3N2)

Note:

Protein SEQ:

qklpgndnstatlclghavpngtiktndrievtnatlvqnssigeicdsphqildgenctlidallgdpqcdgfnkkw
dlfverskaysncypdyvpyaslrslvassgtlefnnesfnwagvtqngtssscirgsnssffsrlnwlthlnskypalnv
pneqfdklyiwgvhphgtdkdqiflyaqssgritvstkrsqavipnigrprirdipsrisiywtivkpgdillinstgnliap
rgyfkirsgkssimrdsapigkcksecitpngsipndkpfqvnritygacpryvkqstlklatgmrnvperqtrg