



## HA1(H3N2)(A/Wuhan/359/95)

---

Product Name	HA1(H3N2)(A/Wuhan/359/95)
Catalog No	IT-003-0049p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(H3N2)(A/Wuhan/359/95)(17-345aa; GenBank# AAT64865)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

