

HA1 (A/Cambodia/R0405050/2007 (H5N1))



Product Name	HA1(A/Cambodia/R0405050/2007(H5N1))
Catalog No	IT-003-00526p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA1(A/Cambodia/R0405050/2007(H5N1)) protein (a.a. 24-341)(Genebank No. ACI06178).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(A/Cambodia/R0405050/2007(H5N1))