

HA(Δ TM)(H5N1)(A/chicken/Viet Nam/NCVD-03/2008)



Product Name	HA(Δ TM)(H5N1)(A/chicken/Viet Nam/NCVD-03/2008)
Catalog No	IT-003-00527 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA(Δ TM)(H5N1)(A/chicken/Viet Nam/NCVD-03/2008) (GenBank No. FJ842481; aa 17-532)
Applications	WB standard, antibody ELISA, etc
Purity	\geq 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H5N1)(A/chicken/Viet Nam/NCVD-03/2008)

Note: