


HA (A/ Bar-headed Goose/ Qinghai/ 12/05(H5N1))



Product Name	HA (A/ Bar-headed Goose/Qinghai/12/05(H5N1))
Catalog No	IT-003-0055p-PBS
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA (A/ Bar-headed Goose/Qinghai/12/05(H5N1)) (Genebank No. ABE68927) protein (a.a. 17-530) with deletion at its cleavage site (Δ RRRKKR). Majority of the protein exists as trimer/oligomer
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	$\geq 95\%$
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified HA (A/ Bar-headed Goose/Qinghai/12/05(H5N1)) protein