




HA₁(H5N1)

Product Name	HA ₁ (H5N1)
Catalog No	IT-003-005p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg purified protein in 100 µl PBS only or PBS with 0.1% BSA and 25% glycerol.
Viral Protein	6xHis tagged Hemagglutinin (amino acid 1-345) (H5N1) (A/Vietnam/1194/2004) (Genebank No. ABP51976)
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



293 cell expressing HA1 (VN1194)

Note: