

HA(Δ TM)(A/chicken/Jalisco/CPA1/2012) (H7N3)



Product Name	HA(Δ TM)(A/chicken/Jalisco/CPA1/2012(H7N3))
Catalog No	IT-003-0071 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(A/chicken/Jalisco/CPA1/2012) (H7N3)(a.a.19-517) (See its amino acid sequence below; Genebank #: AFN85519)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM) (A/chicken/Jalisco/CPA1
/2012(H7N3))

Protein SEQ:

dkiclgghavangtkvntltergievvnatetvetantkkictqgkrptdlgqcglgltligppqcdqflefdadliierregtdvc
ypgkftneeslrqilrgsggidkesmgftysgirtngatsacrrsgssfykemkllnsdnaafpqmtksyrnprnkpalii
wgvhhsgsateqtklygsgnklitvgsskyqqsftpsgarpqvnqsggridfhwllldpndvtftngafiapdrasffrge
sigvqsdvpldsgegdcfhnggtivsslpfqinprtvgkcpvkvqtslllatgmrnvpnpkdrksrhrtrglfgaiagfi
engweglidgwygfrhqnaqgegtadykstqsaidqitgklrnlidktnqqfelidnefseieqqignvinwtrdsmtve
wsynaellvamenqhtidladsemnklyervrkqlrenaedgtgcfeifhkcdqcmesirnntydhqtqyraeslqnr