


HA(Δ TM)(A/environment/California/7451/2010(H7N3))



Product Name	HA(Δ TM)(A/environment/California/7451/2010(H7N3))
Catalog No	IT-003-0072 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA(Δ TM)(A/environment/California/7451/2010(H7N3))(a.a.18-520) (See its amino acid sequence below; Genebank #: AFK75166)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ:

dkiclgghavangtkvntltergievvnatetvetanikkiciqrkrptdlgqcglgtligppqcdqflefdadliierregtdvc
ypgkftneeslrqilrgsggidkepmgftysgirtngatsacrrsgssfyamkwllnsdnaafpqmtksynprnkpalii
wgvhhsgsateqtklygsgnklitvgsskyqqsftpsgtrpqvngqsgridfhwllldpndvtftfngafiapdrasflrge
slgiqsdvpldpgegdcfhsggtiisslpfqinprtvgkcpnyvkqtslllatgmrnvpknktrglfgaiagfiengwegli
dgwygfrhqnaggegaadykstqsaidqitgklnrlidktnqqfelidnefseieqqignvinwtrdsmtewwsynaellv
amenqhtidladsemnklyervrkqlrenaedgtgcfeifhkcdqcmesirnntydhtqyrteslqnriqidpvklssg