


# HA1(H17)(A/little yellow-shoulderedbat/ Guatemala/153/2009(H17N10))



---

Product Name	HA1(H17)(A/little yellow-shoulderedbat/Guatemala/153/2009 (H17N10))
Catalog No	IT-003-022p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(H17) (a.a.18-342) (A/little yellow-shoulderedbat/Guatemala/153/2009(H17N10))(Genebank Accession No. AFC35418.1)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(H17)(A/little yellow-shoulderedbat/  
Guatemala/153/2009(H17N10))