

HA1 (B/Hong Kong/05/1972)



Product Name	HA1 (B/Hong Kong/05/1972)
Catalog No	IT-003-B13p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1 (B/Hong Kong/05/1972)(a.a.16-360) (Genebank #: ABF21280)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

