


HA1(B/Finland/139/2010)



Product Name	HA1(B/Finland/139/2010)
Catalog No	IT-003-B16p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(B/Finland/139/2010) (a.a.18-362) (See below for its aa sequence).
Applications	Protein standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein Sequence:

ICTGITSSNSPHVVKATQGEVNVTVIPLTTTPTKSYFANLKGTKTRGKLCPDCLNCTDL
DVALGRPMC VGTTPSAKASILHEVRPVTSGCFPIMHDRTKIRQLANLLRGYENIRLSTQN
VIDAEKAPGGPYRLGTSGSCP NATSKSGFFATMAWAVPKDNNKNATNPLTVEVPYICAE
GEDQITVWGFHSDNKTPMKNLYGDSNPQKFTSSANGVTTHYVSQIGGFDPQTEDGGLPQ
SGRIVVDYMMQKPGKTGTIVYQRGVLLPQKVWCASGRSKVIKGSPLIGEADCLHEKYG
GLNKSHPYYTGEHAKAIGNCPIWVKTP LKLANGTKYRPPAKLLKERG