

# HA( $\Delta$ TM )(A/California/06/2009) (H<sub>1</sub>N<sub>1</sub>)



---

Product Name	HA( $\Delta$ TM )(A/California/06/2009(H1N1))
Catalog No	IT-003-SW1 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(A/California/06/2009(H <sub>1</sub> N <sub>1</sub> )) protein (a.a. 18-530)(Genebank No. <a href="#">ACP41935</a> ).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

