


# NS3 (HCV/Genotype-1b/Con)



Product Name	NS3 (HCV/Genotype-1b/Con)
Catalog No	IT-004-067Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS. Reconstitute the protein using 100 µl Ultrapure water.
Protein	GST-fused HCV NS3 Genotype-1b consensus immunodominant regions (aa 1193-1460).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



## Protein SEQ:

avdfvpvesmettmrsvftdnssppavpqtqfvaahlhaptgsgkstkvpaaayaaqgykvlvlnpsvaatlfggaymska  
hgvdpnirtgvrnittgapitystygkfladggcsggaydiiicdechstdstilgigtvldqaetagarlvvlatatppgsvtvp  
hpnieevalsntgeipfygkaipietikggrhlifchskkkcdelaaklsglgnavayyrgldvsviptsgdvvvvatdalmt  
gftgdfdsvidentevt