




NP(A/Sydney/5/2021)(H1N1)

Product Name	NP(A/Sydney/5/2021)(H1N1)
Catalog No	IT-003-051Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea.
Protein	C-terminal 6xHis tagged NP(A/Sydney/5/2021)(H1N1) (a.a. 1-498, see below)
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NP(A/Sydney/5/2021)(H1N1)

Note:

Protein AA

MASQGTKRSYEQMETGGERQDTTEIRASVGRMIGGIGRFYIQMCTELKLSDYDGRLIQNSITIERMVLSAFDERRNKYLEEHPSAGKD
PKKTGGPIYRRIDGKWTRELILYDKEEIRRVRWQANNGEDATAGLTHIMIWHSNLNDATYQRTRALVRTGMDPRMCSLMQGSTLPR
RSGAAGAAVKGVGTIAMELIRMIKRGINDRNFWRGENGRTRVAYERMCNILKGGFQTAAQRAMMDQVRESRNPNGAEIEDLIFLA
RSALILRGSVAHKSCLPACVYGLAVASGHDFEREGYSLVGIDPFKLLQNSQVVSMLMRPNENPAHKSQVLVWMAHSAAFEDLRVSSFI
RGKKVIPRGKLSRGVQIASNENVEITMDSNTLELRSRYWAIRTRSGGNTNQKASAGQISVQPTFSVQRNLPFERATIMAAFSGNNEG
RTSDMRTEVIRMMESAKPEDLSFQGRGVFELSDEKATNPVPSFDMNSNEGSYFFGDNAEEYDN