



NP(B/Austria/1359417/2021)

Product Name	NP(B/Austria/1359417/2021)
Catalog No	IT-003-053Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea.
Protein	C-terminal 6xHis tagged NP(B/Austria/1359417/2021) (a.a. 1-560, see below)
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NP(B/Austria/1359417/2021)

Note:

Protein AA

MSNMDIDGMNTGTIDKTPEEITSGTSGATRPIIRPATLAPPSNKRTRNPSPERATTSEDDVGRKAQKKQTPTEIKKSVYNMV
VKLGEFYNQMMVKAGLNDMERNLIQNAHAVERILLAATDDKKTQFQKKKNARDVKEGKEEIDHNKTGGTFYKMRVDD
KTIYFSPIRITFLKEEVKTMKYKTTMGSDGFSLNHMIGHSQMNDVCFQRSKALKRVGLDPSLISTFAGSTVPRRSGATGVAIK
GGGTLVAEAIIRFIGRAMADRGLLRDIKAKTAYEKILLNLKNKCSAPQQKALVDQVIGSRNPGIADIEDLTLARSMVVVRPSV
ASKVVLPISYAKIPQLGFNVEEYSMVGYEAMALYNMATPVSILRMGDDAKDKSQLFFMSCFGAAYEDLRVLSALTGTEFKP
RSALKCKGFHVPAAKEQVEGMGAALMSIKLQFWAPITRSGGNEAGGDGGSGQISCPVFAVERPIALSKQAVRRMLSMNIEG
RDADVKGNLLKMMNDSMAKKTSGNAFIGKMFQISDKNKTNPPIPIKQTPNFFFGRDTAEDYDDLIDY