



NS1 (Wesselsbron virus)

Product Name	NS1 (Wesselsbron virus)
Catalog No	IT-006-0061p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged NS1 (Wesselsbron virus) antigen (a.a. 752-1127) (Genbank Accession No. ABI54474)
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>80% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

